

Using Clinisept+ Skin

Clinisept+ Skin is the UK's leading skin prep and aftercare – all in one.

Clinisept+'s unique and award winning hypochlorous technology delivers unsurpassed levels of skin hygiene in clinic, and reduces client downtimes and enhances results following treatments such as injectables, laser treatments, chemical peels, micro-needling and more.

In Clinic - Apply Clinisept+ onto clean skin to deliver optimal pre-procedure hygiene and antimicrobial protection against complications.

Clinisept+ can be either sprayed directly onto the skin or decanted into a gallipot and then carefully wiped on to the skin using a gauze or pad. Ensure the treatment area and surrounding skin is thoroughly wetted. Clinisept+ can continue to be used throughout procedures to keep the skin clean and soothed. Once the procedure is complete, a wetted gauze can be left to soak on the treated area to calm and cool.

At home - Apply Clinisept+ Skin to minimise discomfort and redness and provide the optimum conditions for the skin to heal.

Spray Clinisept+ Skin directly onto the skin, or spray onto a cotton pad and gently wipe over the area using a circular motion 2-3 times a day or as needed to reduce redness, irritation and keep skin clean and healthy. For maximum care of recovering or irritated skin, saturate a gauze or pad with Clinisept+ and leave it to soak on the affected area for five to ten minutes. For best results, leave skin to dry naturally.

Clinisept+ Skin



Clinisept+ Skin is a high purity hypochlorous solution manufactured to GMP ISO 22716.

Meets BS EN1040, BS EN1276, BS EN13727, BS EN1650, BS EN13624, BS EN14476, BS EN13704.




Clinical Health Technologies Ltd

Hinckley, LE10 3BP
+44(0)1455 247 797

www.cliniseptplus.com

Clinisept+ Skin



**The gamechanger in
the aesthetics and
beauty industry**

Clinisept+ Skin has revolutionised the way we clean and care for the skin in the aesthetics and beauty sector. Its unique formulation delivers maximum protection against complications whilst simultaneously providing the perfect environment for optimised patient outcomes – without causing any irritation or sensitisation.

- Highly effective but gentle antimicrobial
- Dermatologically tested
- Non-irritant, non-sensitising, non-cytotoxic and hypoallergenic
- Non-irritant to eyes
- Skin neutral pH
- Free from alcohol or chlorhexidine

No other product combines similar levels of skin hygiene with complete skin and tissue compatibility. Endorsed by world-leading cosmetic professionals and recommended by the JCCP (Joint Council for Cosmetic Practitioners) and the World Aesthetic Complications Expert Group (ACE).



Clinisept+ Technology

What makes Clinisept+ Skin unique?

Most skin cleansers use harsh chemicals to clean the skin. Even at low concentrations these are still toxic and can actually inhibit skin cell regeneration.

Clinisept+ Skin takes a different approach: it contains a pure high purity and highly stable hypochlorous solution which replicates the chemistry of the human immune system. Clinisept+'s gentle oxidising action provides highly effective microbial protection but with a skin neutral pH and without causing any toxicity or leaving any residue on the skin.

Clinisept+ Skin Soothes and Calms

Because of its complete skin compatibility, Clinisept+ Skin reduces the appearance of redness and swelling and helps to create the optimum environment for healthy skin to recover.



Clinisept+ is used and endorsed by leading practitioners and bodies in the cosmetic and aesthetic sector.

“The application of Clinisept+ prior to any aesthetic treatment is key to enhancing patient safety. Clinisept+ marks a new era in skin preparation and is a vital link in the chain of best clinical practice”.



Helena Collier

MSc BSc PGDp RGN DNC NIP

Author of Aesthetic Complications Expert Group World Infection Control Guidelines

Clinisept+ passes EN1040, EN1276, EN13727, EN13624, EN1650, EN14476 and EN13704.

